

Datasheet: MPC 360 ABW

Nursing Cart

Price ex works Sargenroth:

Net list price [EUR]: 2.476,00

Technical data [mm/kg]:

Height: 1145
Width: 915
Depth: 525
Weight: 43
Usable space of top plate: 550 x 410
Usable shelf space (inside): 470 x 435
Capacity per shelf: 50
Subject to technical changes!



Date: 18.01.2024

Features:

- Cabinet with lamellae door with lock and 2 height adjustable shelves inside
- Waste container, 20 liter (apply knee pressure to open lid)
- Cabinet made of 19,6 mm RESOPAL® HPL plates with duroplastic film
- Stable frame made of RCN proprietary high impact plastic
- The material is scratch resistant, hygienic, UV stable, heat resistant and easy to clean
- Recessed shelf on top of cabinet, removable for easy cleaning
- Push-bar
- Easy to clean
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- 3 year warranty
- 5 year warranty on frame

Color selection cabinet (MPC and MDC series)

