

Datasheet: MPC 372 S

Nursing Cart

Price ex works Sargenroth:

Net list price [EUR]: 2.532,00

Technical data [mm/kg]:

Height:	1470
Width:	1540
Depth:	525
Weight:	49
Usable space of top plate:	550 x 410
Usable shelf space (inside):	525 x 435
Height of top plate:	1090
Capacity per shelf:	50
Subject to technical changes!	



Date: 18.01.2024

Features:

- PLEASE NOTE: figure shows optional bags (WSW), not included in the price!
- Cabinet with 2 swing doors, drawer with lock and height adjustable shelf inside
- Cabinet made of 19,6 mm RESOPAL® HPL plates with duroplastic film
- Stable frame made of RCN proprietary high impact plastic
- The material is scratch resistant, hygienic, UV stable, heat resistant and easy to clean
- Liquid dispenser on slide rail (ESP), standard for half liter bottles (for liter bottles upon request)
- Glove dispenser on slide rail (HSP)
- Recessed shelf on top of cabinet, removable for easy cleaning
- Double hamper
- Low noise and odor containing lid, protected against breaking through overextending
- Lid opening with foot pedal
- Easy to handle docking system to attach hamper to cart
- Cabinet with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- Hamper with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 2 locking
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection cabinet (MPC and MDC series)



Color selection bags



Color selection lids

