

# Datasheet: MPC 470

## Nursing Cart

### Price ex works Sargenroth:

Net list price [EUR]: 1.926,00

### Technical data [mm/kg]:

Height:	1200
Width:	930
Depth:	525
Weight:	59
Usable space of top plate:	810 x 410
Usable shelf space (inside):	780 x 435
Height of top plate:	1090
Capacity per shelf:	45
Subject to technical changes!	



Date: 18.01.2024

### Features:

- Cabinet with 2 swing doors, drawer with lock and height adjustable shelf inside
- Cabinet made of 19,6 mm RESOPAL® HPL plates with duroplastic film
- Stable frame made of RCN proprietary high impact plastic
- The material is scratch resistant, hygienic, UV stable, heat resistant and easy to clean
- Recessed shelf on top of cabinet, removable for easy cleaning
- Push-bar
- Easy to clean
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- 3 year warranty
- 5 year warranty on frame
- PLEASE NOTE: figure shows lockable doors (optional), not included in the price!

### Color selection cabinet (MPC and MDC series)

