

Datasheet: PV10 EM10 Privacy Screen

Price ex works Sargenroth:

Net list price [EUR]: 854,00

Technical data [mm/kg]:

Height: 1800

Width: 2000

Depth of running gear: 505

Subject to technical changes!



Date: 01.01.2024

Features:

- Privacy Screen standard model 1 with extension screen
- Basic model can be connected with extension screen- if required both can also be set at an angle to each other
- Basic model with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Extension screen with 3 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame

Color selection Sure-Chek® fabric

