

Datasheet: RPC 800 Recliner

Price ex works Sargenroth:

Net list price [EUR]: 2.004,00

Technical data [mm/kg]:

Height: 1400
Width: 605
Depth: 1040-1680
Weight: 32
Seat height: 620
Seat width: 490
Seat depth: 455
Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

Features:

- Contoured seat and backrest cushions in blue
- Swing-out and removable armrests
- Height adjustable backrest (7 levels)
- Height adjustable legrest (7 levels)
- 5 single-wheel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- Option for more comfort: Tray that clips on to armrests (RSB)