

**Net price ex works Sargenroth:**

Promotion price; no further discount [EUR]: 3.980,00

**Technical data [mm/kg]:**

Height:	1400
Width:	720
Depth:	1100 - 1680
Weight:	60
Seat height:	640
Seat width:	480
Seat depth:	600
Weight capacity:	120

Subject to technical changes!



Date: 01.01.2024

**Features:**

- Continuously adjustable backrest
- Side parts of the backrest adjustable using grid adjustment (can be folded down for transfer)
- Armrests adjustable using grid adjustment (can be folded down for transfer)
- Continuously adjustable leg rest
- Height-adjustable pillow
- Inclination angle of the seat infinitely adjustable
- Padding made from breathable comfort and AngelFlex gel foam core
- Gel foam core offers optimal balance of gentle pressure relief, support and temperature compensation
- Cover made of coated fabric - wipeable and resistant to disinfectants
- Push-bar
- Anti-tip frame with optimal tilt point, made of high quality stainless steel
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- 3 year warranty
- 5 year warranty on stainless steel frame
- Additional option: serving table (can be pushed on to armrests), price upon request

Color selection coated fabric

