

# Datasheet: SCC 200 FT PP

## Multi-Purpose Shower and Commode Chair

### Price ex works Sargenroth:

Net list price [EUR]: 741,00

### Technical data [mm/kg]:

Height:	1020
Width:	555
Depth:	560
Weight:	14
Seat height:	555
Seat width:	440
Seat depth:	475
Clearance under seat:	455
Lower leg length:	435
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

### Features:

- Soft seat with toilet and hygiene opening made of PUR-foam
- Monobloc backrest with push handles, made of solid plastic
- Slide-out, anti-tip footrest
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

### Color selection PUR-foam



### Color selection plastic backrest / plastic seat

