

Datasheet: SCC 200 S

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 652,00

Technical data [mm/kg]:

Height:	1080
Width:	555
Depth:	565
Weight:	12,5
Seat height:	535
Seat width:	440
Seat depth:	490
Clearance under seat:	420
Lower leg length:	415
Weight capacity:	150

Subject to technical changes!



Date: 01.01.2024

Features:

- Monobloc seat with toilet and hygiene opening, made of solid plastic
- Monobloc backrest with push handles, made of solid plastic
- Slide-out, anti-tip footrest
- Safety-bar
- Drop-down armrest
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

Color selection plastic backrest / plastic seat



cobalt graphite red