

Datasheet: SCC 220 MS

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 869,00

Technical data [mm/kg]:

Height:	1000
Width:	650
Depth:	575
Weight:	12,5
Seat height:	480
Seat width:	540
Seat depth:	490
Clearance under seat:	365
Weight capacity:	175
Subject to technical changes!	



Date: 01.01.2024

Features:

- Monobloc seat with toilet and hygiene opening, made of solid plastic
- Monobloc backrest with push handles, made of solid plastic
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Pail with inside scale
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 175 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

Color selection plastic backrest / plastic seat



cobalt graphite red