

Datasheet: SCC 250

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 539,00

Technical data [mm/kg]:

Height: 995
Width: 555
Depth: 600
Weight: 9,7
Seat height: 500
Seat width: 440
Seat depth: 460
Clearance under seat: 420
Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

Features:

- Ergonomically designed toilet seat, lifts up and can be removed
- Backrest made of "Sunsure" fabric
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

Color selection Sunsure fabric

