

Datasheet: SCC 250 Comfort G

Multi-Purpose Shower-/ Commode Tilt Chair

Price ex works Sargenroth:

Net list price [EUR]: 4.216,00

Technical data [mm/kg]:

Height:	1380
Width:	655
Depth:	1050
Weight:	38
Seat height:	550
Width of the soft seat:	440
Seat depth:	470
Weight capacity:	150
Width between armrests:	525
Lower leg length:	380-490
Subject to technical changes!	



Date: 28.02.2024

Features:

- Ideally suited to assist in the care of people with minimal upper body strength or that are hyperactive
- Infinitely adjustable seating position with gentle movements using a gas pressure spring
- Tilt angle infinitely adjustable from -5° to + 40°
- Soft seat with toilet and hygiene opening, nests on frame and removes easily for cleaning
- PUR-foam padded armrests, upward folding
- Soft backrest made of PUR-foam, removes easily for cleaning
- 3D adjustable headrest
- Individually adjustable, foldable and removable footrests
- Pail with inside scale
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Anti-tip frame with optimal tilt point, made of high quality stainless steel
- Guaranteed NO RUST
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on stainless steel frame