

Datasheet:

SCC 250 Comfort M

Multi-Purpose Shower-/ Commode Tilt Chair

Net price ex works Sargenroth:

Promotion price; no further discount [EUR]: 3.551,00

Technical data [mm/kg]:

Height:	1335
Width:	650
Depth:	1145
Weight:	42
Seat height:	565
Width of the soft seat:	440
Seat depth:	530
Weight capacity:	150
Width between armrests:	515
Lower leg length:	430-505
Subject to technical changes!	



Date: 01.01.2024

Features:

- Ideally suited to assist in the care of people with minimal upper body strength or that are hyperactive
- Hand crank ensures continuous adjustable sitting position and smooth movements
- Tilt angle infinitely adjustable from -5° to + 40°
- Soft seat with toilet and hygiene opening, nests on frame and removes easily for cleaning
- PUR-foam padded armrests, upward folding
- Soft backrest made of PUR-foam, removes easily for cleaning
- 3D adjustable headrest
- Individually adjustable and foldable footrests
- Pail with inside scale
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Anti-tip frame with optimal tilt point, made of high quality stainless steel
- Guaranteed NO RUST
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on stainless steel frame