

Datasheet: SCC 250 D SB

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 909,00

Technical data [mm/kg]:

Height:	1040
Width:	555
Depth:	605
Weight:	16
Seat height (without cushion):	530
Seat width:	440
Seat depth (without cushion):	460
Clearance under seat:	450
Lower leg length:	410
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Safety-bar
- Swing-out armrest (on the right, by request on the left)
- Slide-out, anti-tip footrest
- Ergonomically designed toilet seat, lifts up and can be removed
- Backrest made of "Sunsure" fabric
- Backrest cushion and seat cushion with foam core and cover made of coated fabric
- With cushions also applicable for use as room chair (cushions must to be taken off in shower!)
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- Padding thickness (unloaded): 50 mm (backrest cushion) and 60 mm (seat cushion)
- Compression hardness of the foam: 3 kPa
- 3 year warranty
- 5 year warranty on frame

Color selection coated fabric

