

# Datasheet: SCC 250 FT PR

## Multi-Purpose Shower and Commode Chair

### Price ex works Sargenroth:

Net list price [EUR]: 752,00

### Technical data [mm/kg]:

Height:	1050
Width:	555
Depth:	590
Weight:	15
Seat height:	530
Seat width:	440
Seat depth:	450
Clearance under seat:	450
Lower leg length:	410
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

### Features:

- Ergonomically designed toilet seat, lifts up and can be removed
- Slide-out, anti-tip footrest
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

### Color selection PUR-foam



graphite blue cobalt