

Datasheet: SCC 250 MS

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 887,00

Technical data [mm/kg]:

Height:	985
Width:	660
Depth:	590
Weight:	13
Seat height:	480
Seat width:	540
Seat depth:	470
Clearance under seat:	365
Dimension of toilet opening:	205 x 295
Weight capacity:	175
Subject to technical changes!	




Date: 01.01.2024

Features:

- Monobloc seat with toilet and hygiene opening, made of solid plastic
- Backrest made of "Sunsure" fabric
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 175 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection Sunsure fabric



Color selection plastic backrest / plastic seat

