

Datasheet: SCC 250 OS DS ASAOS PR

XXL Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 1.820,00

Technical data [mm/kg]:

Height:	995
Width:	680-720
Depth:	590
Weight:	19
Seat height:	480
Seat width:	560-600
Seat depth:	450
Clearance under seat:	400
Dimension of toilet opening:	260 x 210
Weight capacity:	250
Subject to technical changes!	




Date: 01.01.2024

Features:

- Reinforced swing-out armrest (on the right, by request on the left)
- Flat plastic seat with toilet opening
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 250 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection PUR-foam



graphite blue cobalt