

Datasheet: SCC 250 OS PPG STE

XXL Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 2.405,00

Technical data [mm/kg]:

Height:	1060
Width:	680-720
Depth:	670
Weight:	31
Seat height:	565
Seat width:	560-600
Seat depth:	540
Clearance under seat:	465
Lower leg length:	400
Dimension of toilet opening:	260 x 210
Weight capacity:	250

Subject to technical changes!



Date: 01.01.2024

Features:

- Seat depth extension to 54 cm
- Reinforced swing-out armrest (on the right, by request on the left)
- Slide-out, roll-along footrest
- Soft seat with toilet and hygiene opening made of PUR-foam, only available in blue
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Footrest with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 250 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame