

Datasheet: SCC 250 OTP SB

Orthopedic Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 2.052,00

Technical data [mm/kg]:

Height:	1040
Width:	555
Depth:	630
Weight:	23
Seat height:	555
Seat width:	440
Seat depth:	450
Clearance under seat:	455
Lower leg length:	435
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Adjustable tilt angle of the backrest and the seat
- Tilt angle of the backrest adjustable from 100° to 120° (with the seat in horizontal position)
- Tilt angle of the seat adjustable from -3° to + 3°
- Soft seat with toilet and hygiene opening made of PUR-foam
- Insert for soft seat
- Two-part soft backrest made of PUR-foam, only in blue
- Safety-bar
- Swing-out armrest (on the right, by request on the left)
- Slide-out, anti-tip footrest
- Pail with inside scale
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- 3 year warranty
- 5 year warranty on frame