

Datasheet: SCC 250 PPG SB AP ER

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 1.195,00

Technical data [mm/kg]:

Height:	1065
Width:	555
Depth:	590
Weight:	21,5
Seat height:	565
Seat width:	430
Seat depth:	450
Clearance under seat:	465
Lower leg length:	405
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Safety-bar
- Swing-out armrest (on the right, by request on the left)
- PUR-foam padded armrests
- Slide-out, anti-tip footrest
- Soft seat with toilet and hygiene opening made of PUR-foam
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 5 year warranty on frame
- 3 year warranty

Color selection PUR-foam



graphite blue cobalt