

Datasheet: SCC 250 PPG SB DRS

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 1.175,00

Technical data [mm/kg]:

Height:	1050
Width:	555
Depth:	590
Weight:	20
Seat height:	555
Seat width:	440
Seat depth:	450
Clearance under seat:	455
Lower leg length:	435
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Safety-bar
- Swing-out armrest (on the right, by request on the left)
- Protection against slipping off the seat (can be easily taken off), PUR-foam padded
- Slipping protection eliminates the need for a safety-belt
- Slide-out, anti-tip footrest
- Soft seat with toilet and hygiene opening made of PUR-foam
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 5 year warranty on frame
- 3 year warranty
- Option for more comfort: insert for soft seat (PSE)

Color selection PUR-foam



graphite blue cobalt