

Datasheet: SCC 250 PPG SB HV

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 1.084,00

Technical data [mm/kg]:

Height:	970-1070
Width:	555
Depth:	585
Weight:	18,4
Seat height:	480-580
Seat width:	440
Seat depth:	450
Clearance under seat:	380-480
Lower leg length:	360-460
Weight capacity:	150
Subject to technical changes!	




Date: 01.01.2024

Features:

- Height adjustment to 5 different levels (25 mm steps)
- Safety-bar
- Swing-out armrest (on the right, by request on the left)
- Slide-out, anti-tip footrest
- Soft seat with toilet and hygiene opening made of PUR-foam
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 5 year warranty on frame
- 3 year warranty

Color selection PUR-foam



graphite blue cobalt