

Datasheet: SCC 250 PP FT

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 806,00

Technical data [mm/kg]:

Height:	1040
Width:	555
Depth:	605
Weight:	15
Seat height:	555
Seat width:	440
Seat depth:	460
Clearance under seat:	455
Lower leg length:	435
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:


- Slide-out, anti-tip footrest
- Soft seat with toilet and hygiene opening made of PUR-foam
- Backrest made of "Sunsure" fabric
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on frame

Color selection PUR-foam



graphite blue cobalt

Color selection Sunsure fabric



white navy blue suncast blue wineberry laguna yellow confetti