

Datasheet: SCC 250 PR PZ

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 783,00

Technical data [mm/kg]:

Height:	1010
Width:	555
Depth:	590
Weight:	11
Seat height:	525
Seat width:	440
Seat depth:	450
Clearance under seat:	425
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Soft seat with continuous hygiene split for total access, made of PUR-foam, color blue
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on frame

Color selection PUR-foam

