

# Datasheet: SCC 250 R

## Multi-Purpose Shower-/ Commode Travel Chair

### Price ex works Sargenroth:

Net list price [EUR]: 1.126,00

### Technical data [mm/kg]:

Height: 815-890  
Width: 555  
Depth: 555  
Weight: 12,3  
Seat height: 470-545  
Seat width: 440  
Seat depth: 470  
Clearance under seat: 390-465  
Weight capacity: 150  
Bag dimension: 630x600x340  
Subject to technical changes!



Date: 06.03.2024

### Features:

- No tools needed to assemble or disassemble
- With travel bag, made of "Sunsure" fabric
- Ergonomically designed toilet seat, lifts up and can be removed
- Backrest, padded with PUR-foam
- PUR-foam padded armrests
- Height adjustment to 4 different levels (25 mm steps)
- Swing-out armrest on both sides
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Pail with inside scale
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- Guaranteed NO RUST
- 3 year warranty
- 5 year warranty on frame

### Color selection PUR-foam



### Color selection Sunsure fabric

