

Datasheet: SCC 250 RC

Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 2.477,00

Technical data [mm/kg]:

Height:	1380
Width:	555
Depth:	1000
Weight:	33
Seat height:	595
Seat width:	440
Seat depth:	430
Weight capacity:	150
Dimensions of backrest (WxH):	420x810
Dimensions of legrest (WxD):	400x400
Subject to technical changes!	



Date: 01.01.2024

Features:

- Height adjustable backrest (7 levels)
- Height adjustable legrest (6 levels) with welded watertight pad
- Adjustable angle between flat seat and backrest from 100° to 146°
- Total length with reclined backrest and legrest in horizontal position increases to 1.630 mm
- Drop-down armrest on the right
- Soft seat with toilet and hygiene opening made of PUR-foam, only available in blue
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 5 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 2 locking
- Pail with inside scale
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 5 year warranty on frame
- 3 year warranty