

# Datasheet: SCC 250 XL PPG STE

## XXL Multi-Purpose Shower and Commode Chair

### Price ex works Sargenroth:

Net list price [EUR]: 1.854,00

### Technical data [mm/kg]:

Height:	1060
Width:	580-660
Depth:	670
Weight:	31
Seat height:	565
Seat width:	460-540
Seat depth:	540
Clearance under seat:	465
Lower leg length:	400
Dimension of toilet opening:	260 x 210
Weight capacity:	200

Subject to technical changes!



Date: 01.01.2024

### Features:

- Seat depth extension to 54 cm
- Swing-out armrest (on the right, by request on the left)
- Slide-out, roll-along footrest
- Soft seat with toilet and hygiene opening made of PUR-foam, only available in blue
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Pail with inside scale
- Lid storage on backrest
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 200 kg
- Easy to clean
- 5 year warranty on frame
- 3 year warranty