

Datasheet: SCC 250 XXXL DSP 300 MFTOS PR STE XXXL Multi-Purpose Shower and Commode Chair

Price ex works Sargenroth:

Net list price [EUR]: 3.915,00

Technical data [mm/kg]:

Height:	1060
Width:	820
Depth:	735
Weight:	38
Seat height:	550
Seat width:	700
Seat depth:	540
Dimension of toilet opening:	260 x 210
Weight capacity:	300
Lower leg length:	385
Subject to technical changes!	




Date: 01.01.2024

Features:

- Slide-out, roll-along footrest
- Flat plastic seat with toilet and hygiene opening
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Seat depth extension to 54 cm
- 6 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 4 locking
- Pail with inside scale
- Lid storage on backrest
- Double framed back made of RCN proprietary high impact plastic, reinforced with inside aluminum tubing
- UV-resistant and compatible with most disinfectants
- Guaranteed NO RUST
- Does NOT roll over toilet
- Supports up to 300 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection PUR-foam



graphite blue cobalt