

Datasheet: SC 100 PPG SB 2ASA HV

Shower Chair

Price ex works Sargenroth:

Net list price [EUR]: 1.192,00

Technical data [mm/kg]:

Height:	975-1075
Width:	555
Depth:	590
Weight:	20
Seat height:	480-580
Seat width:	440
Seat depth:	450
Clearance under seat:	380-480
Lower leg length:	360-460
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2023

Features:

- Swing-out armrest on both sides
- Height adjustment to 5 different levels (25 mm steps)
- Slide-out, anti-tip footrest
- Safety-bar
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Soft seat with toilet and hygiene opening made of PUR-foam
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection PUR-foam



graphite blue cobalt