

Net price ex works Sargenroth:

Promotion price; no further discount [EUR]: 199,00

Technical data [mm/kg]:

Height:	1005
Width:	540
Depth:	580
Weight:	8,5
Seat height:	505
Seat width:	440
Seat depth:	445
Clearance under seat:	440
Weight capacity:	130
Subject to technical changes!	



Features:

- Back cover made of durable Sunsure fabric in blue
- Armrests with soft-grip coating in blue
- Ø 75 mm twin-wheel swivel casters with lock
- Full plastic toilet seat with hygiene opening
- Rolls over toilet
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 130 kg
- Easy to clean
- 2 year warranty

Date: 01.01.2026