

# Datasheet: SDH 70

## Stand Stool

### Price ex works Sargenroth:

Net list price [EUR]: 517,00

### Technical data [mm/kg]:

Height: 910

Width: 450

Depth: 510

Weight: 8,4

Seat height: 650

Seat width: 390

Seat depth: 200

Weight capacity: 120

Subject to technical changes!



Date: 01.01.2024

### Features:

- For people that cannot sit all the way down but need something to rest on
- Seat- and backrest pad made of PUR-foam, only available in blue
- Special anti-slip legs
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 120 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame