

# Datasheet: SM 700 FBS SBS Beach Mobil

## Price ex works Sargenroth:

Net list price [EUR]: 3.539,00

## Technical data [mm/kg]:

Height: 1130

Width: 1065

Depth: 1125

Weight: 40

Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

## Features:

- Rolls easily over the most soft, sandy beaches thanks to its low pressure balloon tires
- Hand operated brake/lock
- Safety-bar
- Swing-out armrest
- Seat- and Backrest Cushion with foam core and cover made of "Sunsure" fabric
- Can be knocked down into three sections, no tools needed
- Tote bag made of "Sunsure" fabric
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Salt water resistant
- Guaranteed NO RUST
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

## Color selection Sunsure fabric

