

Net price ex works Sargenroth:

Promotion price; no further discount [EUR]: 3.814,00

Technical data [mm/kg]:

Height:	1920
Width:	1150
Depth:	635
Weight:	72
Inner dimension cabinet:	490x520x540
Inner dimension drawer (DxWxH):	435x460x110
Inner dimension technique compartment:	555x520x235
Subject to technical changes!	



Date: 18.01.2024

Features:

- Lets you take your senses activation and enjoyment therapy to different places in your facility
- Cabinet with 2 swing doors, 1 adjustable shelf inside, 1 lockable drawer and 1 lockable technique compartment
- Technique compartment opens by lifting the top plate and slide it to the back
- Glass cylinder to hold water for light effects
- Color drapes on rotating and removable rod
- High-quality Hi-Fi system
- Light globe
- Scent atomizer with light effects
- Cabinet made of 19,6 mm RESOPAL® HPL plates with duroplastic film
- The material is scratch resistant, hygienic, UV stable, heat resistant and easy to clean
- Frame made of RCN proprietary high impact plastic
- 4 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Easy to clean
- 3 year warranty (electrical equipment excluded)
- 5 year warranty on frame
- The senses activation cart can also be equipped with technique package 2 (price on request)
- Technique package 2 contains the "Mathmos Space Projector with lava lamp effect" instead of the light globe

Color selection cabinet (MPC and MDC series)

