

Datasheet: TC 450

Transport Chair

Price ex works Sargenroth:

Net list price [EUR]: 1.324,00

Technical data [mm/kg]:

Height: 1125
Width: 555
Depth: 650
Weight: 22
Seat height: 570
Seat width: 440
Seat depth: 450
Weight capacity: 150
Lower leg length: 400

Subject to technical changes!



Date: 01.01.2024

Features:

- Soft seat made of PUR-foam
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Slide-out, anti-tip footrest
- Push-bar
- Swing-out armrest on both sides
- PUR-foam padded armrests
- Safety-bar
- 4 single-wheel swivel casters (plastic/stainless steel), Ø 100 mm, all locking
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame

Color selection PUR-foam

