

Datasheet: TC 450 K

MRI compatible Transport Chair up to 3 Tesla

Price ex works Sargenroth:

Net list price [EUR]: 2.423,00

Technical data [mm/kg]:

Height:	1120
Width:	550
Depth:	650
Weight:	22
Seat height:	565
Seat width:	440
Seat depth:	450
Weight capacity:	150
Lower leg length:	390
Subject to technical changes!	



Date: 01.01.2024

Features:

- Completely made of plastic, without any metal parts
- Suitable for all MRI systems and electric security gates (does not activate an alarm) - MR SAFE
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Soft seat made of PUR-foam
- Slide-out, anti-tip footrest
- Push-bar
- Swing-out armrest on both sides
- PUR-foam padded armrests
- Safety-bar
- 4 single-wheel swivel casters (all-plastic), Ø 125 mm, 2 locking
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 2 year warranty

Color selection PUR-foam



graphite blue cobalt