

# Datasheet: TC 450 MRT

## MRI compatible Transport Chair up to 3 Tesla

### Price ex works Sargenroth:

Net list price [EUR]: 1.453,00

### Technical data [mm/kg]:

Height: 1125  
Width: 550  
Depth: 650  
Weight: 21  
Seat height: 570  
Seat width: 440  
Seat depth: 450  
Weight capacity: 150  
Lower leg length: 420

Subject to technical changes!



Date: 01.01.2024

### Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Soft seat made of PUR-foam
- Soft backrest made of PUR-foam, removable for easy cleaning and access
- Slide-out, anti-tip footrest
- Push-bar
- Swing-out armrest on both sides
- PUR-foam padded armrests
- Safety-bar
- 4 twin-wheel swivel casters (plastic/brass), Ø 100 mm, all locking, non-magnetic
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Supports up to 150 kg
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required

### Color selection PUR-foam

