

Datasheet: TL 500 S

Shower and Transfer Gurney

Price ex works Sargenroth:

Net list price [EUR]: 1.959,00

Technical data [mm/kg]:

Height:	940
Width:	820
Depth:	1970
Weight:	28
Deck height:	765
Up to shoulder width:	595
Up to body length:	1860
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Comfortable sling deck made of robust "Sunsure" fabric
- Drop-down side rails
- Frame made of RCN proprietary high impact plastic
- Easy to clean
- UV-resistant and compatible with most disinfectants
- Guaranteed NO RUST
- Supports up to 150 kg
- 4 single-wheel swivel casters (plastic/stainless steel), all locking
- 3 year warranty
- 5 year warranty on frame

Color selection Sunsure fabric

