

Datasheet: TL 500 S MRT

MRI compatible transfer gurney up to 3 Tesla

Price ex works Sargenroth:

Net list price [EUR]: 2.089,00

Technical data [mm/kg]:

Height:	940
Width:	820
Depth:	1970
Weight:	28
Deck height:	765
Up to shoulder width:	595
Up to body length:	1860
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Comfortable sling deck made of robust "Sunsure" fabric
- Drop-down side rails
- Frame made of RCN proprietary high impact plastic
- Easy to clean
- UV-resistant and compatible with most disinfectants
- Guaranteed NO RUST
- Supports up to 150 kg
- 4 twin-wheel swivel casters (plastic/brass), Ø 100 mm, all locking, non-magnetic
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required

Color selection Sunsure fabric

