

# Datasheet: TL 600 MRT

## MRI compatible transfer gurney up to 3 Tesla

### Price ex works Sargenroth:

Net list price [EUR]: 3.964,00

### Technical data [mm/kg]:

Height:	960
Width:	760
Depth:	2010
Weight:	48
Deck height:	770
Up to shoulder width:	630
Up to body length:	1900
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

### Features:

- Fulfills all safety requirements for stations up to 3 Tesla - MR CONDITIONAL
- Comfortable deck pad with foam core
- Cover made of coated fabric
- Drop-down side rails
- Height adjustable head/backrest (6 levels)
- Frame made of RCN proprietary high impact plastic
- Easy to clean
- UV-resistant and compatible with most disinfectants
- Guaranteed NO RUST
- Supports up to 150 kg
- 4 twin-wheel swivel casters (plastic/brass), Ø 100 mm, all locking, non-magnetic
- 3 year warranty
- 5 year warranty on frame
- For stations higher than 3 Tesla we offer customer specific product solutions, even all plastic equipment if required

### Color selection coated fabric

