

Datasheet: UDL 1500 VR

Shower and Transfer Gurney

Price ex works Sargenroth:

Net list price [EUR]: 3.702,00

Technical data [mm/kg]:

Height:	960
Width:	880
Depth:	2020
Weight:	55
Deck height:	800
Up to shoulder width:	610
Up to body length:	1820
Weight capacity:	150
Subject to technical changes!	



Date: 01.01.2024

Features:

- Height adjustable head/backrest (6 levels)
- Soft tub-liner offers comfort
- Drop-down siderails
- Drain hose, 90 cm
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- Guaranteed NO RUST
- 4 single-wheel swivel casters (plastic/stainless steel), all locking
- Supports up to 150 kg
- 3 year warranty
- 5 year warranty on frame
- Option for more comfort: wedge pillow (DK)