

# Datasheet: UDW 2500 VRW

## Shower and Transfer Gurney

### Price ex works Sargenroth:

Net list price [EUR]: 6.970,00

### Technical data [mm/kg]:

Height: 700-1100

Width: 790

Depth: 2000

Weight: 86

Deck height: 490-890

Up to shoulder width: 620

Up to body length: 1850

Weight capacity: 200

Subject to technical changes!



Date: 01.01.2024

### Features:

- Adjustable backrest
- Soft tub-liner offers comfort
- Drop-down side- and foot-end rails
- Headend rail can be put in a lateral position to extend the length of the deck (extended deck length 2000 mm)
- Height adjustable via footpedals on both sides
- Anti- / Trendelenburg positioning via manually activated pneumatic-spring
- Drain hose, 90 cm
- Easy to clean and to disinfect
- Frame made of powder-coated stainless steel
- Guaranteed NO RUST
- Supports up to 200 kg
- 8 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 4 locking
- 3 year warranty
- 5 year warranty on stainless steel frame
- With aid number in Germany: 04.40.03.4009
- Option for more comfort: wedge pillow (DK)