

# Datasheet: UDW 650

## Shower and Transfer Gurney

### Price ex works Sargenroth:

Net list price [EUR]: 5.620,00

### Technical data [mm/kg]:

Height: 700-1100

Width: 710

Depth: 1770

Weight: 69

Deck height: 490-890

Up to shoulder width: 540

Up to body length: 1620

Weight capacity: 150

Subject to technical changes!



Date: 01.01.2024

### Features:

- Soft tub-liner offers comfort
- Drop-down side- and foot-end rails
- Headend rail can be put in a lateral position to extend the length of the deck (extended deck length 1770 mm)
- Height adjustable via footpedals on both sides
- Anti- / Trendelenburg positioning via manually activated pneumatic-spring
- Drain hose, 90 cm
- Easy to clean and to disinfect
- Frame made of powder-coated stainless steel
- Guaranteed NO RUST
- Supports up to 150 kg
- 8 twin-wheel swivel casters (plastic/stainless steel), Ø 100 mm, 4 locking
- 3 year warranty
- 5 year warranty on stainless steel frame
- Option for more comfort: wedge pillow (DK)