

Datasheet: WK 07 WL

Baby changing unit

Price ex works Sargenroth:

Net list price [EUR]: 8.932,00

Technical data [mm/kg]:

Height:	1770
Width:	1010
Depth:	860
Weight:	92
Top of changing pad:	970
Surface of changing pad (WxD):	825 x 760
Inner dimension drawer (WxDxH):	240 x 580 x 160
Length of bathtub (inside):	580
Width of bathtub (inside):	260
Height of bathtub (inside):	200
Subject to technical changes!	



Date: 01.01.2024

Features:

- Padded changing surface with foam core and cover made of coated fabric
- Changing table with three side padded rail (200 mm height from top of the changing surface)
- High quality upholstery is compatible with most disinfectants and resistant to urine and blood
- Heat lamp holder including heat lamp unit (type: IW 980 AGU, Fisher & Paykel)
- Heat lamp unit with skin temperature control
- Heat lamp can be swung to the side
- Distance between changing surface and lower edge of the heat lamp: 680 mm
- Integrated cabinet with 1 storage compartment and 2 drawers
- Cabinet made of 19,6 mm RESOPAL® HPL plates with duroplastic film
- Bathtub cart with 4 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, all locking (height of cart: 875 mm)
- Bathtub color blue, including drain hose (hose length 840 mm)
- Changing unit with 6 twin-wheel swivel casters (plastic/stainless steel), Ø 75 mm, 4 locking
- Frame made of RCN proprietary high impact plastic
- UV-resistant and compatible with most disinfectants
- Easy to clean
- 3 year warranty
- 5 year warranty on frame
- 1 year warranty on heat lamp

Color selection bathtub



light blue gray

Color selection cabinet (MPC and MDC series)



leaf-green gentian blue signal red clementine orange broom yellow alabaster white

Color selection coated fabric



saffron cherry sky pine charcoal